



Nitri-Care® # 3005PF

100% NITRILE GLOVE
 DISPOSABLE 4 MIL, 9.5" GLOVE
 CLASS ONE MEDICAL GLOVE
 US PATENT 35.616



Disposable Glove:

- Powder-Free
- Silicone-Free
- Non-Sterile Exam Gloves
- Made in China
- Food Handling: In Compliance 21 CFR, Parts 170 to 199

Description

Best's Nitri-Care® 3005 Class One Medical Examination glove is a value-based option for those seeking disposable hand protection solutions in medical and Emergency Medical Transport applications as well as food packaging and handling. Powder-free and medically approved as a Class One Medical Device, Nitri-Care® comes in a pebbled finished to ensure excellent grip and handling in a wide range of sizes.

Suggested Applications

- Patient Examination
- Manufacturing
- Food Processing
- Laboratories
- Janitorial
- Paint Applications
- Fuel Handling
- Venipuncture

This list of suggested applications is partial and based on laboratory testing of glove performance. For a complete list of suggested applications, visit our website:

www.showabestglove.com

Part Number	Size	Packaging
3005PFXS	Extra-Small	100 Gloves per Dispenser 10 Dispensers per Case
3005PFS	Small	
3005PFM	Medium	
3005PFL	Large	
3005PFXL	Extra-Large	



USA:
 Showa Best Glove
 579 Edison Street
 Menlo, GA 30731
 Phone: 800-241-0323
 Fax: 888-393-2666
 Email: USA@showabestglove.com

Canada:
 Showa Best Glove Manufacturing, Ltd.
 253 Michaud Street
 Coaticook, Quebec J1A 1A9 Canada
 Phone: 800-565-BEST
 Fax: 819-849-6120
 Email: Canada@showabestglove.com

Europe:
 Showa Best Glove
 Tour Franklin
 La Defense 8
 92042 Paris La Defense Cedex France
 Tel Direct +33 (0)1 55 62 14 20
 Fax + 33 (0)1 55 62 14 29
 Email: Info@showabestglove.eu

Australasia:
 Best Manufacturing Australia Pty. Ltd.
 PO Box 187, Bexley
 NSW 2207 Australia
 Phone: +61 (0)2 8509 5401
 Fax: +61 (0)2 8509 5403
 Email: Australiansales@bestglove.com