

SemperSure™

NITRILE EXAMINATION GLOVES

Industrial

Exam

Surgical

Non-Latex • Powder-Free • Textured



Sempermed® 

We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



SemperSure™ by Sempermed nitrile examination gloves are designed for those with skin sensitivities to regular nitrile gloves. SemperSure™ accelerator free gloves offer all the benefits of a nitrile glove, minus the sulfur-based chemical accelerators that cause Type I (latex) and Type IV (chemical) allergies.

FEATURES

- Accelerator free
- Textured surface on fingers
- Tested and approved for use with chemotherapy drugs
- Manufactured without natural rubber
- Powder-free

BENEFITS

- Low dermatitis potential
- Ensures a safe grip, even in a wet environment
- Safe for use with chemotherapy drugs
- Eliminates exposure to latex
- Eliminates aerosolization of glove powder

GLOVE COUNT (By Weight)

200 Gloves/Box (180 Gloves/Box XL)
10 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.05 mm/2.0 mils minimum
Palm: 0.05 mm/2.0 mils minimum

PHYSICAL PROPERTIES*

Before Aging:

Tensile Strength: 14 MPa min.

Elongation: 500% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.

All Sempermed gloves are single use only.

PRODUCT CHARACTERISTICS

- Synthetic Nitrile Rubber (NBR)
- Accelerator Free
- Low Dermatitis Potential
- Tested for Use with Chemotherapy Drugs
- Textured
- Powder-Free
- Ambidextrous
- Single Use Only

Order #	Size	Length (min.)	Palm Width (min.)
SUNF201	XS	240mm	75mm +/- 5mm
SUNF202	S	240mm	85mm +/- 5mm
SUNF203	M	240mm	95mm +/- 5mm
SUNF204	L	240mm	106mm +/- 5mm
SUNF205	XL	240mm	116mm +/- 5mm

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558